

Extremely Low Magnification EBSD with an FE-SEM

Y. Yamamoto*, K. Ogura*, K. Kawauchi*, M. Shibata*, N. Erdman**, and C. Nielsen**



* JEOL Ltd., 3-1-2 Musashino, Akishima, Tokyo, Japan

** JEOL USA, Inc. 11 Dearborn Rd, Peabody, MA 01960, U.S.A.


Electron Backscattered Diffraction (EBSD) with an SEM has recently become very popular for the crystal orientation analysis of materials. It is quite useful to analyze tiny crystals at rather high magnification from several hundreds up to several tens thousands times. On the other hand, for the analysis of very large crystals with a size of several hundreds micrometers, a very low magnification SEM/EBSD is required. However, the very low magnification SEM/EBSD is somewhat problematic, because of the distortion of SEM images due to a large tilting angle of specimens. Recently, we have developed an FE-SEM (JSM-7001F/TTL), which resolved the image distortion problem at low magnification. SEM/EBSD with this JSM-7001F/TTL is significantly less affected by the image distortion even in the magnification range of less than one hundred times.

JSM-7001F/TTL equipped with an EBSD system (TSL OIMTM) was used to investigate a sample of polysilicon for solar cell application. The JSM-7001F/TTL is a high-resolution Schottky FEG-SEM. For high-magnification imaging, an electron probe is focused to a specimen by the objective lens (OL). On the other hand, when the low magnification imaging is necessary, the electron probe is focused to the specimen using the Aperture Angle Control Lens (ACL) instead of the OL. In this Low Magnification (LM) mode, the working distance between the ACL to the specimen is much longer as compared to using the OL. Furthermore, the aperture angle of the electron probe becomes much smaller as well. As the result, the LM mode of the JSM-7001F/TTL has an enormous depth of field without image distortion, thus allowing low-magnification SEM/EBSD with minimal artifacts. Fig. 1 shows a comparison between the conventional imaging mode and the LM mode. Fig. 2 shows an image comparison at ten times of the magnification between the conventional imaging mode and the LM mode using a 70 degree tilted metal grid.

Fig.3 shows an optical microscope image of solar cell polysilicon with the size of 30 mm square. The area marked with a black rectangle line was used for SEM/EBSD. The actual SEM magnification for SEM/EBSD was x12. The analysis was performed at 20 kV, 5 nA of electron probe current, working distance 25 mm, 70 degree of specimen tilt, and 100 micron EBSD measurement step. Fig. 4 shows results of low magnification EBSD, such as an image quality (IQ) map, and an inverted pole figure (IPF) maps. These results demonstrate that the LM mode of the JSM-7001F/TTL is suitable for the extremely low magnification SEM/EBSD applications.



Normal mode



Low Mag. mode

Fig.2: SEM image of a Cu grid

Fig.1: Electron optics of JSM-7001F/TTL




Fig.3: Polysilicon for solar cell




Fig. 4: An image quality map and inverse pole figure maps