

# An E.T. Tube That Talks Back To You.

## **Smooth Murphy Eye**

Provides an alternate ventilation port should the tip of the tube become occluded.

## **High Volume, Low Pressure Blue Cuff.**

Blue cuff communicates a visual landmark for tube placement and provides an adequate cuff seal to the tracheal wall at minimum occluding pressures. This reduces the risk of restricting blood flow in the tracheal wall, minimizing the risk of necrosis or tissue breakdown.

## **Gently Cupped Tip**

The distal portion of the MatrX tube has a gently cupped tip which provides an easy, non-traumatic transition of the tube past the vocal chords.

## **Bold, Easy To Read Depth Marks**

Visually tells you the proper tube placement depth in low light situations.

## **Soft Yet Solid**

Though soft, the tube provides adequate strength to accomplish intubation.

## **Precision Engineered Inflation System**

Includes the highly reliable Bepak inflation valve. The "pilot balloon" reacts quickly and provides excellent communication between you and the cuff. You can immediately feel pressure changes within the cuff after intubation. This will allow you to determine, by the touch of your hand, the relative pressure the cuff is exerting on the tracheal wall. This minimizes the potential for tracheal tissue damage. Upon extubation, the "pilot balloon" will deflate to a completely flat position when the air is withdrawn from the cuff.

**Matrx**  
MEDICAL INC.

*Quality and pride in all we do.*

MATRIX IS A FULL LINE DEALER. U.S. ORDERS CALL 800-845-3550

CANADA 800-847-1000 • INTERNATIONAL ORDERS 716-662-6650 • INTERNATIONAL FAX: 716-662-7130

# FROM RESUSCI ANNE TO THE HEART OF THE MATTER

