

INDEX TO ADVERTISERS

Diatome	Cover 3
EDAX	A7
Electron Microscopy Sciences, EMS	Cover 2
Hitachi	Cover 4
Physik Instrumente	A4
Tescan	A8

The best way to get better results



The APEX™ platform is the powerful and user-friendly software that drives the collection and analysis of Energy Dispersive Spectroscopy (EDS) and Electron Backscatter Diffraction (EBSD) data for the compositional and microstructural characterization of materials.

Please visit edax.com/apex for more information.

Faster results from an intuitive user interface

Customizable color schemes and layouts

Best analytical performance for both EDS and EBSD

Easy reporting for the flexible presentation of data


AMETEK®

MATERIALS ANALYSIS DIVISION

edax.com

EDAX®
Smart Insight

Take gallium-free microsample preparation further



Array of micro-compression test pillars in UFG aluminium alloy

TESCAN AMBER X

- ✓ High throughput, large area FIB milling up to 1 mm
- ✓ Ga-free microsample preparation
- ✓ Ultra-high resolution, field-free FE-SEM imaging and analysis
- ✓ In-column SE and BSE detection
- ✓ Spot optimization for high-throughput, multi-modal FIB-SEM tomography
- ✓ Superior field of view for easy navigation
- ✓ Essence™ easy-to-use, modular graphical user interface



For more information visit


www.tescan.com

TESCAN AMBER X is based on the S8000 platform.

DiATOME Diamond Knives



at the forefront of innovation



NEW! Ultra ATS The Ultra ATS Diamond knife is perfect for placing sections on Si wafers to view under the SEM. The knife comes in 3.0mm size with 35° angle.

**ultra 45° • cryo • histo • ultra 35°
histo jumbo • STATIC LINE II • cryo immuno
ultra sonic • ultra AFM & cryo AFM
trimtool 20, 45, and 90**

www.emsdiasum.com

DiATOME U.S.

P.O. Box 550 • 1560 Industry Rd. • Hatfield, PA 19440
Tel: (215) 412-8390 • Fax: (215) 412-8450

email: info@emsdiasum.com or stacie@ems-secure.com

Size without Compromise

The SU3900 VP-SEM with an Extra-Large Chamber

Hitachi's SU3900 Variable-Pressure SEM features an extra-large chamber to accommodate a multitude of analytical applications: EDS, EBSD, micro XRF, in situ, heating, cooling, tensile testing, and more. Equipped with 20 ports and a robust stage, the SU3900 can handle it all! Whether your sample is large, small, heavy, magnetic, or non-conductive, the SU3900 is the instrument for your lab.

Extra-Large Chamber & Heavy-Sample Stage

Superior X/Y range of motion for up to 5 kg sample weight

Ultra-Wide NaviCam FOV

Covering full observation sample area (200 mm Φ)

Intelligent Filament Technology

Auto optimization for improved filament performance and lifetime

Fully Integrated High-Quality EDS Capability


Seamless SEM and EDS operation from a single GUI



For more information:
microscopy@hitachi-hta.com



Sample: Pure Titanium
Accelerating Voltage: 5 kV
Magnification: 800x



Extra-large chamber for max. 300 mm ϕ , 130 mm tall, and 5 kg sample

Innovation • Synergy • Solutions



Science for a better tomorrow

©Hitachi High-Tech America, Inc.

www.hitachi-hightech.com/us Tel. 800-253-3053 E-mail: microscopy@hitachi-hta.com