


Microscopy TODAY



Volume 27 Number 6 2019 November



The Hitachi ArBlade 5000 Advanced Ion-Milling System

The Hitachi ArBlade 5000 Ion-Milling System is a highly advanced broad ion-beam system. The ArBlade 5000 equipped with a fast-milling Ar ion gun will increase your throughput and put your lab on the cutting edge.

**PUT YOUR LAB
ON THE CUTTING EDGE**



- 1 **New-generation hybrid instrument with cross-section milling and Flatmilling™**
- 2 **≥ 1,000 μm per hour cross-section milling on Si—Cut your milling time in half!**
- 3 **Cross-section widths up to 8 mm—Mill out massive areas with the new cross-section holder**

Visit Hitachi at 2019 MRS Fall Meeting
Dec. 1-6, 2019, Boston, MA

Innovation ♦ Synergy ♦ Solutions

Expand your Knowledge of Microscopy with MSA Membership!

Whether your primary focus is in light, electron, or scanning probe microscopy, or the biological or the physical sciences, MSA takes your knowledge to the next level!


Members Receive:

- A personal subscription to MSA's official journal, *Microscopy and Microanalysis*, and MSA's popular bi-monthly magazine, *Microscopy Today*.
- Peer Networking through the Society's Focused Interest Groups and Local Affiliated Societies.
- Discounts on books, journals and other educational materials.
- MSA Awards Programs, Scholarships, Speaker Opportunities, and more!



Join MSA Today!
Visit www.microscopy.org

First-and-Only AFM-Based Viscoelastic Technique that Ties Directly to Bulk DMA



AFM-nDMA mode – Addressing Nanoscale Domains with Highest Resolution

How do nanoscale phases differ from bulk? What are viscoelastic properties right near an interface or in an interphase? How does this affect load transfer, as a function of temperature? Bruker’s new AFM-nDMA mode allows researchers to answer these types of questions. For the first time an AFM can provide complete and quantitative viscoelastic analysis of polymers at the nanoscale, probing materials at rheologically relevant frequencies, in the linear regime. Proprietary dual-channel detection, phase-drift correction, and reference frequency tracking enable a small strain measurement in the rheologically relevant 0.1 Hz to 20 kHz range for nanoscale measurements of storage modulus, loss modulus, and loss tangent that tie directly to bulk DMA.

NEW AFM-nDMA mode

- Provides accurate moduli and loss tangent at bulk DMA frequencies – a first for AFM
- Delivers absolute calibration with quantified load and adhesion – and no need for a reference sample
- Creates master curves with WLF and Arrhenius analyses from AFM data
- Discovers nanomechanical effects that conventional AFM approaches might miss

Your research deserves the very best that AFM technology has to offer!


Visit www.bruker.com/AFM-nDMA, email productinfo@bruker.com or call +1.805.967-2700 for more information.

Contents

Scanning Probe Microscopy

- 12** The ORNL Lectures on Scanning Probe Microscopy, Part 1: Piezoresponse Force Microscopy and Spectroscopy of Ferroelectrics, Energy Materials, and Biological Systems
Sergei V. Kalinin
- 18** The ORNL Lectures on Scanning Probe Microscopy, Part 2: The Force Dimension: Electronic and Ionic Transport Measurements via Kelvin Probe Force Microscopy
Sergei V. Kalinin
- 24** Correlative Microscopy and Nanofabrication with AFM Integrated with SEM
Mathias Holz, Christoph Reuter, Ahmad Ahmad, Alexander Reum, Martin Hofmann, Tzvetan Ivanov, and Ivo W. Rangelow
- 32** Microscopy 101: Scanning Probes or Scanning Electrons: A Practical Guide to Select a Method for Nanoscale Characterization
Dalia Yablon and Matthew Libera

About the Cover



Multivariate analysis of tubular spinel-perovskite interfaces. Maps of individual transport components are clockwise from upper left: linear Ohmic component, non-linear component, low-conductivity component, and memristive component.

See article by Kalinin, Part 2.

Specimen Preparation

- 40** Preparing TEM Specimens and Atom Probe Tips by Laser Machining
Boris Rottwinkel, André Kreuzer, Henry Spott, Michael Krause, Georg Schusser, and Thomas Höche

Microscopy Pioneers

- 46** Pioneers in Optics: Albert A. Michelson
Eric Clark

Highlights from *Microscopy and Microanalysis*


- 52** Exploring the Parameter Space of Point Spread Function Determination for the Scanning Electron Microscope – Part I: Effect on the Point Spread Function
- 52** Annexin A2 Enhances the Progression of Colorectal Cancer and Hepatocarcinoma via Cytoskeleton Structural Rearrangements
- 53** Electron-excited X-ray Microanalysis by Energy Dispersive Spectrometry at 50: Analytical Accuracy, Precision, Trace Sensitivity, and Quantitative Compositional Mapping

Departments


- | | |
|--------------------------------------|--------------------------------|
| 7 Editorial | 54 NetNotes |
| 8 Carmichael's Concise Review | 62 Calendar of Meetings |
| 48 Industry News | 68 Dear Abbe |
| 50 Product News | 70 Index of Advertisers |

Providing Microscopy Supplies and Specimen Preparation Equipment to Our Valued Customers for Over Half a Century


Visit us December 1–6 in Boston, MA
MRS FALL | BOOTH 616



Adhesives
Conductive Tape, Carbon Tabs
Carbon, Silver & Gold



Calibration Standards,
Pelcotec™ CDMS-XY, Magnification,
Resolution, LM/SEM/FESEM



NEW Quantifoil®
TEM Substrates



**NEW Products for Specimen
Thinning in Preparation for
Electron Microscopy**



Tools for Microscopy
Tweezers, Scribes, Vacuum Pick-Up



**PELCO® Modular SEM
Holders & Mounts**



FIB Supplies
Lift Out Grids, Storage,
Probes, FIB Mounts



Cressington Coating Systems
Carbon Evaporation & Metal Sputtering



**NEW PELCO® Precision
Low Speed Saw**
Designed to cut a wide variety
of materials with minimal
subsurface damage



PELCO® AFM Tools & Supplies
AFM Discs, Disc Carriers, AFM Probes,
Diskpenser, Disc Grippers & Pickup Tools



www.tedpella.com sales@tedpella.com 800-237-3526